

Urocap IV



Simply flows

...Innovating the way Urodynamics is delivered!



UROCAP IV: Uroflowmetry Refined

Urocap IV is the latest in a line of reliable, standalone, cost effective and easy to use Uroflowmetry solutions from LABORIE.

As a wireless, battery-powered, water-protected device, it is capable of providing patient privacy and functional versatility; simply place it where you want and it will send data wirelessly to your computer.

The UROCAP IV software is unparalleled in its ease of use; combined with i-LIST™ Reporter software, you can review and generate a variety of scripted or data reports with a few clicks of a mouse! With one extra click, the reports can be printed, sent to MS Word, PDF writer or automatically emailed to referring physicians. Our latest technology allows you to simplify your clinical workflow and increase patient throughput without strain on your team or compromise on quality of care.

As your solution provider for quality Uroflowmetry, LABORIE's Urocap IV is designed with your work flow in mind.



GO with LABORIE

Technical Specifications

Size & Weight

- L = 6.3" (16.0cm)
- W = 6.3" (16.0cm)
- H = 2.6" (6.6cm)
- Weight = 0.7lbs (0.3kg)

Measurement Channels

- Volume
 - Range = 0 to 1200ml
 - Accuracy = +/- 2%
 - Sampling Rate = 20Hz
- Flow
 - Range = 0 to 50mL/s
 - Accuracy = +/- 5%
 - Sampling Rate = 20Hz

Indicator Lights

- Bluetooth, Status & Battery

Charger

- 5 VDC, 3.0A

Optional Accessories



Commode Chair



Uroflow Stand

USA:

Tel.: +1 802 857 1300
800 522 6743
Fax: +1 802 878 1122

CANADA:

Tel.: +1 450 671 5901
888 522 6743
Fax: +1 450 671 7182

UK:

Tel.: +44 1275 372332
Fax: +44 1275 375145

FRANCE:

Tel.: +33 1 43 41 34 22
Fax: +33 1 43 41 38 69

INTERNATIONAL:

Tel.: +1 905 612 1170
Fax: +1 905 612 9731

BRAZIL:

Tel.: +55 11 4134 8480
Fax: +55 11 4208 3837

BELGIUM:

Tel.: +32 3 780 17 38
Fax: +32 3 780 17 39

CHINA:

Tel.: +86 10 8489 8891
Fax: +86 10 8489 8511



www.laborie.com